

**USP <797> Chapter 2.3 Competency Testing in Aseptic Manipulation**

For additional Cultivate™ Directions for Use and test logs go to [www.parasolmed.com/cultivate](http://www.parasolmed.com/cultivate)

**Definitions**

**Negative:** Light to medium amber, clear, no precipitates

**Positive:** Characterized by turbidity, haziness, specks, or flocculation in the medium.

**TV100 - TSB Growth Media**

**NEGATIVE**



**POSITIVE**



*Escherichia coli*




*Staphylococcus aureus*

**TV020 - TSB Growth Media**

**NEGATIVE**



**POSITIVE**



*Escherichia coli*



*Staphylococcus aureus*